



Data Management System (DMS) *Pro* User Guide

Eversense and the Eversense logo are trademarks of Senseonics, Incorporated. Other brands and their products are trademarks or registered trademarks of their respective holders.

Eversense DMS does not contain malware or computer viruses that might harm your computer. It is recommended you use an appropriate firewall program and an anti-virus program on your computer to prevent against malware and viruses.

To read the Eversense Privacy Policy, go to global.eversenseddiabetes.com.

Table of Contents

Glossary.....	1	4. DMS Pro Dashboard.....	16
I. Introduction.....	2	Viewing Patient Eversense CGM Data	16
Intended Use	2	Setting Date Range	17
End User Licence Agreement and Privacy Policy	2	AGP in DMS Pro	18
Help and Support	2	5. Reports.....	21
2. Getting Started.....	3	Time in Target	23
Eversense DMS Pro Administrator	3	Time in Range	24
Registering a clinic and creating your admin account	3	Glucose Variability Report	25
Logging in to your Eversense DMS Pro Account	4	Glucose Trend Report	26
Eversense DMS Pro Clinic Staff User	8	Glucose History Report	27
Your Profile and the Navigation Menu	8	Glucose Distribution Report	28
Profile	9	Daily Glucose Report	29
Forgot Your Password	10	Transmitter Log Report	31
Home Page	10	6. Settings.....	32
3. Permissions.....	11	Glucose	32
View	12	7. Help.....	33
Add	13	About	33
Edit	15	8. Troubleshooting.....	34
Admin	15		

Glossary

AutoSync When this feature is turned on in the Eversense mobile app, your glucose data and settings will automatically get synced to the Eversense DMS every two hours.

Blood Glucose The glucose value obtained from a fingerprick test.

Blood Glucose Meter A commercially available device used to measure glucose using a blood sample from a fingerprick.

Calibration Blood glucose reading from a fingerprick sample entered in the Eversense App to check the accuracy of the system.

CGM Continuous Glucose Monitoring. Continuously monitoring your glucose levels from interstitial fluid every few minutes.

Contraindication A condition or circumstance in which a person should not use the device.

EULA End User Licence Agreement.

Eversense App Software program that is installed on a mobile device and displays CGM glucose data obtained from the smart transmitter.

Eversense Client Application A desktop application that allows Eversense users to upload the data on their smart transmitter to the Data Management System using a USB cable.

Eversense Data Management System (DMS) A web-based application compatible with the Eversense App where your glucose data is stored and can be viewed.

Interstitial Fluid (ISF) The fluid between cells in the body. The Eversense CGM measures glucose from an interstitial fluid sample, versus glucose in a blood sample obtained from a fingerprick.

mg/dL Milligrams per decilitre, a unit of measure that shows the concentration of a substance in a specific amount of fluid. In some countries, including the United States, glucose test results are reported as mg/dL, indicating how much glucose is in the blood when using a blood glucose meter, or how much glucose is in the interstitial fluid when using some CGM systems, like the Eversense CGM System.

mmol/L Millimoles per litre, a unit of measure that shows the concentration of a substance in a specific amount of fluid. In many countries, glucose test results are reported as mmol/L, indicating how much glucose is in the blood when using a blood glucose meter, or how much glucose is in the interstitial fluid when using some CGM systems, like the Eversense CGM System.

Sensor Glucose The glucose value obtained from the Eversense Sensor.

Sensor A device inserted subcutaneously for continually measuring interstitial fluid glucose levels.

Smart Transmitter A reusable device worn externally over the inserted sensor that powers the sensor and sends glucose information to the mobile device for display in the Eversense App.

I. Introduction

Intended Use

The Eversense Data Management System (DMS) Pro is a tool for healthcare professionals to manage and review glucose data of their patients that use the Eversense CGM System.

For information on compatible browsers and browser versions, please go to <https://global.eversenseddiabetes.com/compatibility>.

End User Licence Agreement and Privacy Policy

Use of the Eversense DMS Pro is subject to the terms and conditions of the most current Eversense End User Licence Agreement and Eversense Privacy Policy. These documents are updated from time to time and are posted at global.eversenseddiabetes.com.

Help and Support

For questions regarding content in this DMS Pro User Guide, contact Customer Support toll-free in the US at 844-SENSE4U (844-736-7348). Outside the US, call your local distributor or visit global.eversenseddiabetes.com to locate your local distributor.

To access the latest version of the Eversense DMS PRO Instructions for use, please visit: global.eversenseddiabetes.com/hcpuserguides.

2. Getting Started

2 Eversense DMS Pro Administrator

Eversense DMS Pro requires at least one account administrator who can:


- Register the clinic, which generates the unique clinic ID.
- Invite other staff members to join the clinic site and set up their permissions.
- Change staff member permissions.
- Delete staff members from the clinic admin account.
- Assist in password reset for other staff members.
- Enable and disable multi-factor authentication and mandatory password expiration for the clinic.

Registering a clinic and creating your admin account

To register your clinic, you must also register a clinic admin account. To do so, go to globalpro.eversensedms.com. After registering, your clinic will be assigned a Clinic Number which will be used as your identifier for patients to be able to request to join your clinic.

To register your clinic and create your admin account:

1. Go to globalpro.eversensedms.com
2. Click on “Sign Up” on the login screen



3. Fill out registration form

IMPORTANT: You cannot use the same email address as a username for Eversense DMS Pro and a personal Eversense DMS account.

4. Click “Submit”
5. You will receive an email from Eversense DMS when your clinic and admin accounts are approved. The email will contain a link to activate your account.

The link will expire after 72 hours.

6. Once you click the link to activate your account, you will be prompted to create your Eversense DMS Pro password

After your password is created, you can log in to your account.


Logging in to your Eversense DMS Pro Account

To access Eversense DMS Pro, go to globalpro.eversensedms.com. You will be directed to the login page after selecting your region.


2

1. Go to globalpro.eversensedms.com
2. On the login page, select your language, enter your account email and password, and choose “Sign In”. You will be directed to your DMS Pro account Home Page.


IMPORTANT: If you make more than three unsuccessful attempts to sign in, your Eversense DMS Pro account will be locked for 10 minutes.



The screenshot shows the login page for healthcare professionals. At the top, it says "For Healthcare Professionals" in a purple banner. Below that is the Eversense logo and "Data Management System Pro". There is a language selector for "English UK". The form includes fields for "Email" (with "sample@email.com" entered) and "Password". A "Sign In" button is highlighted with a red box. Below the button are links for "Need an Eversense DMS Pro Account? Sign Up", "Forgot your Password?", and "Problem logging in? Contact Us". At the bottom, it says "To sign in to Eversense DMS for your CGM System click here" and the Senseonics logo.

3. At the top left-hand side of the Home page is the Navigation Menu. You can open and close the Navigation Menu by clicking . You will also see your Clinic Number below your clinic name.

Note: You can also select your preferred language by clicking the flag icon in the top right corner of the page.




To add staff members:

1. Click “Staff” on the left Navigation Menu
2. Click “Add Staff”
3. Fill out their contact information, and select their title
4. By default, the staff member title determines the permissions assigned. Review and edit the permissions as needed, according to the descriptions below.

Multiple permissions may be assigned.

- **Admin:** Can View, Add, and Edit.
- **View:** Can view the accepted patients list and view the data in each patient’s DMS account.
- **Add:** Same permissions as View, plus the ability to accept patient requests to join your clinic, and to send patients invitations to join your clinic.
- **Edit:** Can edit a patient’s information in DMS Pro, unsubscribe a patient from the clinic, and view their Eversense DMS data.



5. Click “Send”

6. The clinic staff member will receive an email with an account verification link and temporary password. After successfully logging in, they will be prompted to create a new password.

Note: If a staff member cannot find the invitation email you sent when adding them as staff, you have the option to resend the invitation by clicking the “Resend Invitation” button on the **Edit clinic staff** screen (see *next page*). It is recommended to have them check their spam or junk folder.

To edit staff member information and permissions:

LAST NAME	FIRST NAME	TITLE	EMAIL	STATUS	EDIT	RESET PASSWORD
lawrence	john	Physician Assistant	johnlaw12@mailinator.com	Active	Edit	Reset

1. Click “Staff” on the left Navigation Menu
2. Click “Edit” to open their information page
3. From this page you can:
 - a. Edit their name and/or title
 - b. Update their permissions
4. Click “Update” when finished

IMPORTANT: If the title of a user is changed, you should check the permissions and edit as needed.

To deactivate clinic staff accounts:

1. Click “Staff” on the left Navigation Menu
2. Click “Edit” to open their information page
3. Choose Inactive on the bottom of the page
When an account is inactive, the clinic staff will no longer be able to add, edit, or view any patient accounts in DMS Pro.
4. Click “Update” on the bottom right of the page

Edit clinic staff ✕

Contact Information

First Name *

Last Name *

Title *

Active Inactive

Permission *

✓	Permission
✓	Admin
✓	Add
✓	Edit
✓	View

Resend Invitation Update Cancel

Note: To filter the staff list by account status, use the filter option on the top right of the Clinic Staff Accounts page.

Clinic Staff Accounts

Search Filter by:

LAST NAME	FIRST NAME	TITLE		
lawrence	john	Physician Assistant	john	om

To assist in password reset for clinic staff:

LAST NAME	FIRST NAME	TITLE	EMAIL	STATUS	EDIT	RESET PASSWORD
lawrence	john	Physician Assistant	johnlaw12@mailinator.com	Active	Edit	Reset

1. Click “Staff” on the left Navigation Menu
2. In the **RESET PASSWORD** column, select “Reset” across from the name of the appropriate staff member
3. Enter the new password and confirm it
4. Click notify user if you wish to notify the staff member via email that their password has been reset
5. Click “Reset Password”

Reset clinic staff password ✕

Username

New password

- Must contain 6-32 characters
- At least one lower case character
- At least one upper case character
- At least one number
- At least one special character

Confirm Password

Notify user

Reset Password Cancel

To enable multi-factor authentication and/or mandatory password expiration:

Multi-factor authentication helps protect your clinic accounts by requiring a level of identity verification beyond basic login credentials. The default is for multi-factor authentication to be disabled, but the DMS Pro Clinic Administrator can enable multi-factor authentication for all clinic accounts. If enabled, every 365 days all staff members who are Eversense DMS Pro users will need to verify their identity prior to logging into their individual account. They will do this using a verification code sent to the email used to create their account.


Mandatory password expiration helps protect your clinic accounts by requiring users to routinely update their account passwords. Based upon the period of time you select, all staff members within your Eversense DMS Pro clinic will be prompted to create a new password to log into their individual account every 365 days or every 180 days.

1. On the top right of the Home page, click your user name

A drop-down will appear displaying your profile picture, your most recent log in, and options to view/edit your profile or sign out.

2. Click “Profile” to access User Profile Settings


3. Click “Account Settings”



4. Navigate to the bottom of the Profile Details section and check the box next to “Check to enable multi-factor authentication” and/or “Check to enable mandatory password expiration”

The image shows the 'Profile Details' form. It includes fields for Profile Picture (Choose File), First Name (Natalie), Last Name (Mellon), Mobile Number, Date Of Birth, Secret Question (What is the name of your first pet?), and Answer. At the bottom of the form, there are two checkboxes: 'Check to enable multi-factor authentication' (unchecked) and 'Check to enable mandatory password expiry' (checked). Below the second checkbox are two radio buttons: 'Every 365 Days' (selected) and 'Every 180 Days'. The entire checkbox area is highlighted with a red box.

5. On the corresponding pop-up, choose “Enable”



6. At the bottom of the page, click “Save” to save your changes. If you enabled multi-factor authentication, the next time a clinic staff member logs into their account, they will be prompted to complete the multi-factor authentication steps. If you enabled mandatory password expiration, clinic staff members will be prompted to reset their password in either 365 days or 180 days, based upon the time period you selected.

The User Profile Settings page will appear again. To return to the Home page, click “Home” at the top left on the Navigation Menu.

To disable multi-factor authentication and/or mandatory password expiration:

1. Navigate to the bottom of the Profile Details section and uncheck the box next to “Check to enable multi-factor authentication” and/or “Check to enable mandatory password expiration”

2. At the bottom of the page, click “Save” to save your changes

The User Profile Settings page will appear again. To return to the Home page, click “Home” at the top left on the Navigation Menu.

Eversense DMS Pro Clinic Staff User


2

Your admin will create an account for you in Eversense DMS Pro. You will receive an email with a temporary password and a link to activate your account and to create a permanent password. Click on the link to complete your account set up, and display the Home page. Once your account is activated, you can go to globalpro.eversensedms.com to log in at any time.

Your Profile and the Navigation Menu

Upon successful sign in, the Home page will be displayed.

- On the top right of the page you can access your **Profile** information.
- On the top left side of the page you will see a **Navigation Menu**. Your Navigation Menu may look different depending on the permission assigned to you by your clinic admin.



Profile

Located at the top right of the Home page, this feature allows you to change your name or password, edit your user profile, or sign out of the Eversense DMS Pro.

2

Change your password:

1. Click your username

A drop-down will appear displaying your profile picture, your most recent log in, and options to view/edit your profile or sign out.



2. Click “Profile” to access User Profile Settings

3. Click “Change Password” to display and change your current password

4. Enter your Current Password

5. Enter and confirm your new password

6. Click “Change Password” to save your changes

A screenshot of the 'Change Password' form. The form has three input fields: 'Current Password', 'New password', and 'Confirm new password'. The 'New password' field has a list of requirements: 'Must contain 6-32 characters', 'At least one lower case character', 'At least one upper case character', 'At least one number', and 'At least one special character'. At the bottom right of the form, there are two buttons: 'Back to Profile' and 'Change Password' (highlighted with a red box).

To adjust your account information:

1. Click “Account Settings”

2. Enter the information you will be changing

You can change your profile picture by pressing the **Choose File** button, and choosing a picture from your computer's files.



3. Click “Save” to save your changes

The User Profile Settings page will appear again. To return to the Home page, click “Home” at the top left on the navigation panel. Patients using Eversense DMS or clinic staff users using Eversense DMS Pro will not be able to view your profile picture.

A screenshot of the 'Profile Details' form. It features a profile picture placeholder with a 'Choose File' button (highlighted with a red box) and the text 'no file chosen'. Below this are several input fields: 'First Name' (filled with 'Natalie'), 'Last Name' (filled with 'Mellon'), 'Mobile Number', 'Date Of Birth', 'Secret Question' (filled with 'What is the name of your first pet?'), and 'Answer' (filled with '.....'). At the bottom, there are two checkboxes: 'Show Answer' (unchecked), 'Check to enable multi-factor authentication' (unchecked), and 'Check to enable mandatory password expiry' (unchecked). Under the last checkbox, there are two radio buttons: 'Every 365 Days' (selected) and 'Every 180 Days'.

Forgot Your Password

1. Click “Forgot your password” on the Sign In page
2. Enter your email when prompted
3. Click “Send Mail”






A message will display to let you know that an email has been sent to you with instructions for resetting your password.

4. Open the email and click the website link

You will be directed to a web page to create and confirm a new password.

5. Click “Reset Password”

A message will appear, confirming your password has been reset. Select the “click here” link to return to the Sign In page.

Home Page

The default setting for the Home page is to display the list of active patients. Until patients have been added to your DMS Pro clinic, that page will remain blank. Refer to the *Add (Data Sharing)* section that follows. Once populated, you can sort the data by clicking on any column header.

- The List of Patients can be filtered by the provider using the “**Patients by Provider**” button on the top right of the screen.
- Using the drop-down, you can also select to view the inactive patient list. (For more on the inactive patient list, see *section 3, To remove a patient from the clinic*).




PATIENT LAST NAME	PATIENT FIRST NAME	MRN	GE	PROVIDER NAME	LAST SYNC DATE	STATUS	EDIT	DELETE PATIENT
Garcia	Jean	12121	Female	Linda Miller	25-Aug-2025	Active	Edit	Delete
Philip	William		Female	Linda Miller	NA	Active	Edit	Delete

3. Permissions


3

In Eversense DMS Pro, there are different permissions a clinic staff user can have within a clinic. The permissions are **Admin**, **View**, **Edit**, and **Add**. The DMS Pro admin of your clinic assigns the permissions to you. The Navigation Menu may look different depending on the permission level assigned to you by the clinic admin.



Search, Filter, Export Patient List

All permission levels are able to search, filter, and export patient lists.



Click			
To	Search	Filter	Export
How	You can search for patients by name in the search bar above the List of Patients.	You can also filter the List of patients by their healthcare provider by clicking Patients by Provider. <ol style="list-style-type: none"> 1. Click "Patients by Provider" 2. Check the box of the provider whose patients you would like to view 3. Click "Proceed" 	You can export the patient list as an Excel sheet when you click "Export Data" .

View

To view a patient's data:

1. Click on the patient's name

3

Mellon Endocrine
Clinic ID: 7536

List of Patients Patient was successfully removed from the clinic

Filter by: Active

PATIENT LAST NAME	PATIENT FIRST NAME	MRN	GENDER	DOB	PROVIDER NAME	LAST SYNC DATE	STATUS	EDIT	DELETE PATIENT
Garcia	Jean	12121	Female	01-Jun-2007	Linda Miller	25-Aug-2025	Active	Edit	Delete
Philip	William		Female	04-Aug-2006	Linda Miller	NA	Active	Edit	Delete

2. The patient's DMS information will be displayed

For more information, see *Viewing Patient Eversense CGM Data*.

Patient Name: Jean Garcia Age: 23 Gender: Female MRN: 12121 jeangarc97@gmail.com Frank Carl, NP

Last 7 Days Last 14 Days Last 30 Days Last 60 Days Last 90 Days Custom Range

Eversense Wear Time : 100%



Glucose Metric	Value
Average Glucose Goal: <154 mg/dL	139 mg/dL
Glucose Management Indicator (GMI) Goal: <7%	6.6%
Glucose Variability Defined as percent coefficient of variation Goal: <36%	29.2%

Add

To add an Eversense CGM user to your patient list, you can either send them an invitation through DMS Pro, or they can request to join your clinic. The patient must already have an Eversense DMS account. Only Eversense CGM users can be added to your Eversense DMS Pro patient list. For more information on how a patient can create a personal Eversense DMS account, they can visit global.eversensedidiabetes.com and view the Eversense DMS User Guide.


Once the patient has created a personal account, you can send them an invitation, or they can send a request to join your clinic.

To send a patient invitation:

1. Click Patients 
2. Click Invite Patients 
3. Type in the email address that the user has used to register their Eversense account

Note: An invitation cannot be sent if the email has not been registered with Eversense DMS.

4. Click “Search Eversense” to confirm if that email address has been registered with Eversense




The screenshot shows a web form titled "Invite patient". It has a search bar for "Email" with the text "sample@email.com" entered. Below the search bar is a blue button labeled "Search Eversense" which is highlighted with a red rectangular box. There is also a "Reset" button below it.

5. If the email has been registered, you will be prompted to enter patient information

- a. To enter a MRN (optional)
- b. To select a provider within your clinic for this patient
- c. Edit their first or last name
- d. Enter their date of birth (optional)
- e. Enter their gender (optional)

6. Once complete, click “Send”

An email invitation will be sent to the patient.




The screenshot shows the "Invite patient" form with several fields filled out: "Email" is "jeangarc97@gmail.com", "MRN" is "12121", "First Name" is "Jean", and "Date Of Birth" is "07-oct-2000". A dropdown menu for "Provider" is open, showing a list of names: "john lawrence", "Natalie Mellon", "Maria Moore", and "Linda Miller". At the bottom left, a blue button with a checkmark and the text "Send" is highlighted with a red rectangular box. A "Back" button is also visible.

7. You will be asked if you want to send another invitation. If you click No, you will be returned to the patient list screen, where you will see the patient you just invited

8. Until the patient has accepted the invitation, Pending Patient Consent will be displayed in the STATUS column

9. Once they accept the invitation, their status will change to Active and the LAST SYNC DATE column will display the date for their most recent glucose data



The screenshot shows a web interface for "Mellon Endocrine Clinic ID: 7536". It displays a "List of Patients" table. The table has columns for Patient Last Name, Patient First Name, MRN, Gender, DOB, Provider Name, Last Sync Date, Status, Edit, and Delete Patient. Two rows are visible: one for "Garcia, Jean" with status "Active" and "Last Sync Date" "25-Aug-2025", and one for "Haney, Patricia" with status "Pending Patient Consent" and "Last Sync Date" "NA". The "Last Sync Date" and "Status" columns for the Garcia patient are highlighted with red boxes.

PATIENT LAST NAME	PATIENT FIRST NAME	MRN	GENDER	DOB	PROVIDER NAME	LAST SYNC DATE	STATUS	EDIT	DELETE PATIENT
Garcia	Jean	12121	Female	01-Jun-2007	Linda Miller	25-Aug-2025	Active	Edit	Delete
Haney	Patricia			19-Aug-1998	Natalie Mellon	NA	Pending Patient Consent	Edit	Delete

Accepting a patient request

Patients can send a request to join your clinic from their personal Eversense DMS account or directly from their mobile app if they have the Eversense 365 CGM System. You will need to provide them with your Eversense DMS Pro clinic number. To view a list of pending requests, click .


3

To accept a patient request:

1. Click Patients


2. Click Pending Requests

You will see a number next to Pending Requests. This number shows you how many requests are currently in the list. If there is no number, then there are no pending requests.



List of Pending Patients							
Search Pending Patient							Export Data
PATIENT LAST NAME	PATIENT FIRST NAME	MRN	GENDER	DOB	PROVIDER NAME	STATUS	DATA SHARING
Philip	William		Female	04-Aug-2006	Natalie Mellon	Pending	Accept Decline

3. Click “Accept” to add them to your clinic



- When you click “**Accept**”, you will be prompted to assign them to a provider from your clinic.
- If you click “**Decline**”, their request will be removed, and they will not be added to the clinic list.


The Eversense CGM user will be notified via email when you accept or reject their request.


Edit

The **Edit** permission allows you to view and edit a patient's information in DMS Pro, unsubscribe a patient from the Eversense DMS Pro clinic list. Any changes made in DMS Pro will not change any information the patient has saved in their personal account.

3


To edit a patient's information:

1. Click Patients 
2. In the **EDIT** column, select **"Edit"** across from the name of the appropriate patient
3. This will open up a window with their information
You can edit all of their information, except for their email address.
4. Click **"Update"**



To remove a patient from the clinic:


1. Click Patients 
2. In the **DELETE PATIENT** column, select **"Delete"** across from the name of the appropriate patient.



PATIENT LAST NAME	PATIENT FIRST NAME	MRN	GENDER	DOB	PROVIDER NAME	LAST SYNC DATE	STATUS	EDIT	DELETE PATIENT
Garcia	Jean	12121	Female	01-Jun-2007	Linda Miller	25-Aug-2025	Active	Edit	Delete
Philip	William		Female	04-Aug-2006	Linda Miller	NA	Active	Edit	Delete

3. On the pop-up, click **"Yes"** to continue
4. The patient will now be moved to the inactive patient list and their Eversense data will no longer be available to view. If you wish to invite the patient to share data again, choose the **"Invite Patient"** button under the **INVITE TO SHARE AGAIN** column.

With the Edit permission, you are also able to accept Pending Requests from patients who want to join the clinic. See *To accept a patient request* for more information.



PATIENT LAST NAME	PATIENT FIRST NAME	MRN	GENDER	DOB	PROVIDER NAME	STOPPED SHARING	INVITE TO SHARE AGAIN
Gracia	Jean			14-Dec-2005	Natalie Mellon, Physician Assistant	01-Jun-2024	Invite Patient

Note: If a patient chooses to discontinue data sharing with your clinic, they will also be moved to the inactive patient list. You will receive an email informing you that they are no longer sharing their data.

Admin


The **Admin** permission allows you all of the same permissions as **View**, **Add**, and **Edit**, plus you are able to add, edit, and remove staff members. To learn more about the **Admin** permissions, See *Eversense DMS Pro Administrator*.

4. DMS Pro Dashboard

Viewing Patient Eversense CGM Data

The Eversense DMS Pro Dashboard is the default page that you will see each time you select a patient in DMS Pro. This page is very similar to what the patient sees when they log in to their personal Eversense DMS account.

4



At the top of the Home page their name, age, gender, MRN, email and assigned health-care provider is displayed.

You can exit this page at any time by clicking the  in the top right corner of the window or Close at the bottom right of the window.


The Home page contains several sections:

- Date Range Selector
- Short cuts to print or email reports (see section 5)
- Glucose Settings Selector (see section 6)
- Eversense Wear Time
- Ambulatory Glucose Profile (AGP)
- Glucose Reports
- System Information, Date Sensor Linked and Days Since Sensor Linked
- Patient Glucose Settings

Until the patient has uploaded sensor glucose data from the Eversense mobile app or their smart transmitter, no data will be displayed.

IMPORTANT: If the patient has autosync turned on, their data will be automatically uploaded to DMS from their mobile app. They can also manually upload their glucose data from their smart transmitter via a USB cable. For more information on data upload, visit eversenseddiabetes.com to view the Eversense DMS User Guide.

IMPORTANT: Any changes made to the settings or date ranges in DMS Pro WILL NOT affect the settings of the patients, Eversense DMS account or Eversense mobile app.



Setting Date Range

At the top of the Home page, you can view reports based on a predefined date range, or a custom date range you set.

To adjust the date range settings:

1. If you would like to view a set of data that is not within one of the predefined date ranges, click the button labelled “Custom Range” to display a drop down calendar.

4




This will display a pop-up that allows you to select your Custom Date Range.

Note: Data accessible through DMS Pro goes as far back as one year

2. Select the desired report *Start Date*

3. Select the desired report *End Date*

4. Click “Ok” when complete



Eversense Wear Time

Eversense Wear Time is the percentage of time the smart transmitter was worn, over the date range the patient selected. 100% is 23.5 hours of daily smart transmitter wear time. The percentage does not include times in which the glucose is not displayed. For example: the beginning of initialisation phase or when an alert is activated that prevents glucose from being displayed.

Eversense Wear Time : 100%

System Information: Located at the bottom of the AGP Dashboard is the System Information section.

This section shows the patient's:

- Transmitter Serial Number
- Sensor Serial Number
- Date sensor linked and Days since linked
- Glucose settings: Glucose settings can be changed in the DMS as well as through the Eversense mobile app (see *Settings* for more information).

System Information	
Transmitter : 701066	
Sensor Serial Number : 643C2A873C1D3C17	
Date Sensor Linked : 09-Nov-2025	
Days Since Sensor Linked : 63	
Glucose Settings	
Low Target : 85 mg/dL	Low Alert : 70 mg/dL
High Target : 180 mg/dL	High Alert : 200 mg/dL

AGP in DMS Pro

AGP in Eversense DMS Pro

The Eversense DMS Dashboard's primary feature is the AGP (Ambulatory Glucose Profile) report, which displays glucose data in three sections.

4

Patient Name: Jean Garcia Age: 24 Gender: Female MRN: 12121 jeangarc97@gmail.com Natalie Mellon

Last 7 Days Last 14 Days Last 30 Days Last 60 Days Last 90 Days Custom Range

[Glucose Settings](#)

Eversense Wear Time : 100%

Jean Garcia DOB: 07-Oct-2000

14 Days: 01-Feb-2025 - 14-Feb-2025

Time In Range

Goals for Type 1 and Type 2 Diabetes

Very High	1%	12%	Goal: <5%
High	11%		Goal: <25%
Target	86%		Goal: >70%
Low	2%	2%	Goal: <4%
Very Low	<1%		Goal: <1%

Each 1% time in range = about 15 minutes

Glucose Metrics

Average Glucose	138 mg/dL
Goal: <154 mg/dL	
Glucose Management Indicator (GMI)	6.6%
Goal: <7%	
Glucose Variability	29.2%
Defined as percent coefficient of variation	
Goal: <36%	

Ambulatory Glucose Profile

AGP is a summary of glucose values from the report period, with median (50%) and other percentiles shown as if they occurred in a single day.

Daily Glucose Profile

Each daily profile represents a midnight-to-midnight period.

Section 1: Glucose Statistics and Targets

Glucose Ranges and Targets: Displays the patient's calculated values within each glucose range as a percentage of device readings over a 24-hour day. Percentages can be converted to minutes/hours in that each 1% time in range equals about 15 minutes.


- **Time in Target:** Ideally, glucose values are in target range $\geq 70\%$
- **Low/Very Low:** Goals – low $< 4\%$; very low $< 1\%$
- **High/Very High:** Goals – high $< 25\%$; very high $< 5\%$

4

Eversense Wear Time: The percentage of time the smart transmitter was worn over the date range selected. 100% is 23.5 hours of daily smart transmitter wear time.

Number of Days in Report and Dates: The AGP defaults to the last 14 days, but a custom date range can be selected (see *Setting Date Range* earlier in this section).

- **Average Glucose (mean):** All glucose values added together, divided by number of readings
- **GMI (glucose management indicator):** Calculated from average glucose; approximates future lab A1c
- **Glucose Variability:** How far apart (wide) the highest and lowest glucose values are; ideally a low number less than or equal to 36% [percent coefficient of variation]



Source: <http://www.agpreport.org/agp/agpreports>


Section 2: Ambulatory Glucose Profile: All glucose profiles from the selected date range are combined to make a one day (24-hour) picture. Ideally, lines stay within the **green outlined area (target range)**.

Heavy line: Median (middle) line where half of the glucose values are above and half are below; ideally, the line is mostly flat and inside the green target range area of 70-180 mg/dL (3.9-10.0 mmol/L)


Darker shading: Area shaded darker shows 50% of the glucose values; ideally, space between is narrow

Lighter shading: 5% of values are above (95% top line) and 5% are below (5% bottom line); ideally, the closer the light shaded area is to the darker shaded area, the better

4



Section 3: Daily Glucose Profiles: Each box is a single day's glucose pattern. They will only display if the date range is set at 14 days or less.



5. Reports


Reports provide various ways to organise and visualise glucose readings over time.

Note: If you select a report and no data is displayed, check that you have selected a date range in which data is available.

IMPORTANT: If autosync is not turned on, data is not being automatically sent to DMS from the Eversense mobile app.

5 Print or Email what is displayed in the reports

On the top right of the Home page are options for printing and emailing the Home page reports.



Click		
To	Email Report	Save Report as PDF
How	Follow pop-up window prompts for emailing a printer-friendly report to a designated recipient.	Click the icon to convert the report into a PDF that allows you to print the report.


Report Preferences

When exporting or printing reports, the default setting is for you to receive all available reports for the last 14 days of data. However, you may customise your preferences.

To customise your preferences:

1. From the top left of the patient Dashboard, select your preferred date range tab
2. Then, on the top right, select either the Email icon () or Export to PDF icon (), to open the **Report Preferences Screen**.
3. Click to select (or de-select) reports
 - For the Time in Range report, the desired time increments must also be selected.
 - For the Daily Glucose Report, the date range cannot exceed 30 days for printing. You can also select to include event icons and the event table.
4. When you have made your selections, you have the option to **“Export to PDF”** or **“Export & Save as Preferences”**. If choose **“Export & Save as Preferences”**, these selections will be saved and become the new default report preferences for your account.


Note: You can also set/change your report preferences in your Account Profile – see *section 2, Profile*.



Report Options



To view the different report options, click on the name of the report in the Report Panel which is located below the AGP report on the Dashboard.

5

When selected, some reports offer additional filtering. By clicking  next to the report name you can select event types days of the week and other options. This filter is not available for the AGP, Time in Range by Time of Day, or Time in Target reports.

Note: If you select a report and no data is displayed, check the following:

- Expand the date range you have selected. If there is no data over the entire date range over the past 90 days, it may be that the patient does not have autosync turned on in their Eversense CGM mobile app. If autosync is not turned on, ask the patient to turn on autosync.
- Ask the patient to also do a manual upload by clicking “**Start Sync**” in their Eversense CGM mobile app.



Time in Target

The Time in Target report shows a pie chart that analyses the percentage of time your patient's glucose is below, within, or above the target range the patient sets in the Glucose Settings.


Viewing the report:

The green wedge represents the percentage of time the patient spent within the target range.

The grey wedge represents the percentage of time the patient spent above the target range.

5

The red wedge represents the percentage of time the patient spent below the target range.



Time in Range

Time in Range by Time of Day

The Time in Range by Time of Day report shows a set of charts that analyse the percentage of time your patient's glucose is within preset glucose ranges.


- Very Low: <54 mg/dL or <3.0 mmol/L
- Low: 54-70 mg/dL or 3.0-3.9 mmol/L
- Within Range: 70-180 mg/dL or 3.9-10.0 mmol/L
- High: 180-250 mg/dL or 10.0-13.9 mmol/L
- Very High: >250 mg/dL or >13.9 mmol/L

5

There are two sections of this report:

Time in Range


This displays your patient's overall percentage of time spent in the preset glucose ranges, based on the date range you have selected.



Time in Range by Time of Day

This section of the report breaks out the overall Time in Range chart into various segmented times of day. You can view your patient's data in segments of 4 hours, 6 hours, or 8 hours. This is adjusted by selecting an option in the top right of the report.

Each segment shows a percentage of time your patient's glucose was within those same segments of Very Low, Low, Within Range, High, and Very High.




Glucose Variability Report

The Glucose Variability Report displays information in two ways:

The Trend Graph:

- Displays summary glucose readings over a 24-hour period.
- Includes multiple days of data overlaid on the same graph.
- Shows the maximum, minimum, and average glucose reading for every hour.
- Provides a visual snapshot of variations in your patient's daily glucose levels.

Note: Place your cursor above a particular average sensor reading to view details about readings during the hour.



From 06:00 AM to 07:00 AM
 Average: 170 mg/dL
 Min: 117 mg/dL
 Max: 255 mg/dL
 10th Percentile: 125 mg/dL
 90th Percentile: 235 mg/dL

The Table:

- Organises glucose data by 4 hour time period.
- Provides a detailed statistical analysis for each period.
- Includes the percentage of readings that are within, above, and below the glucose target range.
- Provides a breakdown of insulin usage and carbs entered for those same time periods.

Note: If you are having trouble displaying the report, try limiting your date range. Unusually large amounts of data may not allow the report to function.

GLUCOSE	00:00-04:00	04:00-08:00	08:00-12:00	12:00-16:00	16:00-20:00	20:00-00:00	TOTAL
% Below Low Alert	1%	2%	1%	2%	2%	1%	1%
% Above High Alert	3%	-	1%	2%	2%	-	1%
Number of Values Below Low Alert	8	10	6	11	16	5	56
Number of Values Above High Alert	21	-	8	12	15	-	56
% Within Target	81%	86%	81%	89%	88%	90%	86%
% Below Target	2%	2%	1%	2%	5%	1%	2%
% Above Target	17%	12%	18%	9%	8%	9%	12%
Average Glucose	149 mg/dL	135 mg/dL	147 mg/dL	130 mg/dL	138 mg/dL	134 mg/dL	139 mg/dL
Lowest Glucose	46 mg/dL	57 mg/dL	60 mg/dL	53 mg/dL	55 mg/dL	48 mg/dL	46 mg/dL
Highest Glucose	342 mg/dL	226 mg/dL	261 mg/dL	346 mg/dL	350 mg/dL	240 mg/dL	350 mg/dL
Standard Deviation	46.4 mg/dL	35.8 mg/dL	40.9 mg/dL	41.9 mg/dL	39.8 mg/dL	32.9 mg/dL	40.4 mg/dL
Number of Values in the 10-90 Percentile	504	527	539	545	540	545	3,189
Total Number of Glucose Values	624	654	672	672	672	672	3,966
Number of Days with at least One Value	13	14	14	14	14	14	14
Average Number of values per Day	48	47	48	48	48	48	283

Glucose Trend Report

The Glucose Trend Report:


- Displays individual glucose readings over a 24-hour period.
- Displays multiple days of data overlaid on the same graph.
- Provides an easy way to identify variability among the days of the week.

Viewing the report

5

Data points in the report are colour-coded with each day of the week displayed in a different colour. You can use the report filter to choose the day/s you want to display on the graph.

Note: Place your cursor above a sensor glucose reading to view details about that reading.



Glucose History Report


The Glucose History Report displays historical CGM and manually entered data in the form of a trend graph. The Glucose History Report provides an efficient way to track sensor glucose readings and calibrations.

Viewing the report

Each data point is a sensor glucose reading. The events that were logged by the patient in the mobile app are displayed in the report at the time they were logged. If you select a date range of more than 30 days, events such as meals, insulin, etc. are not displayed.

5

Note: Place your cursor above a sensor glucose reading to view details about that reading. You may also place your cursor above an event icon to view details about that event.




Glucose Distribution Report

The Glucose Distribution Report:

- Organises your patient’s glucose data according to their glucose targets and alert ranges and displays the data in pie chart format.
- Displays separate pie charts for the current and previous periods.
- Target ranges are set under Glucose Settings.
- Provides a visual check for the number of readings within target during the predefined range selected and a change over time.

5 Viewing the report

The Glucose Distribution Report will display two pie charts based on the predefined date range you select. It will show the current and previous set of data for the date range selected. This data allows you to see a comparison of current to past changes in your patient’s management.



Daily Glucose Report


The daily glucose report displays glucose readings over a 24-hour period for multiple days in a sequence.

The Daily Glucose Report:

- Has a fixed glucose target range of 70-180 mg/dL (3.9-10.0 mmol/L).
- Shows high and low glucose alert levels (red lines) as set by the CGM user.
- Shows key statistics for each day based upon glucose target and alert settings.
- Includes an option to display event icons on the trend graph, and to expand a table of events for each day.

Viewing the report

The Daily Glucose Report displays multiple days in sequence for the date range selected, and provides a way to identify glucose patterns across multiple days. Statistics for each day appear to the right of the trend graph and are based on the standard glucose target range of 70-180 mg/dL (3.9-10.0 mmol/L) and the high and low glucose alerts levels set by the patient. To review additional CGM information for each day (alerts, calibrations, events entered into the mobile app), you can turn on the display of event icons on the trend graph and expand the event table beneath the trend graph to display more information. When expanded, the event table lists information in 4 hour increments aligned with the 4 hour segments in the trend graph.



Statistics:

The statistics table provides analysis of the percentage of time your patient's glucose is within preset glucose ranges.

- Very Low: <54 mg/dL or <3.0 mmol/L
- Low: 54-70 mg/dL or 3.0-3.9 mmol/L
- Within Range: 70-180 mg/dL or 3.9-10.0 mmol/L
- High: 180-250 mg/dL or 10.0-13.9 mmol/L
- Very High: >250 mg/dL or >13.9 mmol/L

Daily average glucose, standard deviation and number of high and low alerts are also included.


Icons and Events:

The icons and events are displayed based upon your patient's glucose alert settings, calibrations, and any events entered into their mobile app.

Note: If desired, you can adjust the CGM alert settings for the report in DMS Pro. (The glucose target settings are fixed at 70-180 mg/dL (3.9-10.0 mmol/L) and cannot be changed.) Adjusting the glucose alert settings will change the placement of the red alert lines on the daily CGM trend graphs, but not the daily statistics. This may assist in decisions regarding the optimal setting for a high or low glucose alert.



Note: When multiple events are recorded in a single cell in the events table, hovering over the cell will display more information.



Transmitter Log Report

The Transmitter Log Report provides a comprehensive data review over the selected date range.

The Transmitter Log Report displays three separate tables:

- **Blood Glucose:** This displays a list view of all manually entered blood glucose and related events in chronological order.
- **Event Entries:** This displays patient events entered through the Eversense Mobile App.
- **Sensor Glucose:** This displays glucose values picked up by the sensor.

5

Blood Glucose					
Search					
Manual Used Not Used					
Date	Time	Blood Glucose	Unit	Calibration	
07-May-2024	06:18 PM	100	mg/dL	Used	
07-May-2024	06:18 PM	100	mg/dL	Used	
07-May-2024	06:11 PM	100	mg/dL	Used	
07-May-2024	06:11 PM	100	mg/dL	Used	
07-May-2024	05:13 PM	70	mg/dL	Not Used	
07-May-2024	05:13 PM	70	mg/dL	Not Used	
07-May-2024	01:05 PM	100	mg/dL	Manual	
07-May-2024	11:55 AM	93	mg/dL	Used	
07-May-2024	11:53 AM	100	mg/dL	Manual	
07-May-2024	11:53 AM	93	mg/dL	Manual	

Viewing the report

Events are displayed with information regarding Date, Time, and:

- **Blood Glucose:** Glucose Value, Unit, Calibration
- **Event Entries:** Event Type, Event Sub Type, Value
- **Sensor Glucose:** Glucose Value, Unit, Sensor ID

Events can be sorted in ascending or descending order by clicking the column header. The number of records is indicated at the bottom of each log.

A symbol under the “Calibration” column of the Blood Glucose table will indicate whether a glucose value was

 Used for calibration,  Not Used, or  Manual.

- **Used:** A blood glucose value that your patient entered in the Eversense mobile app as a calibration and is accepted by the system.
- **Not used:** A calibration that was entered that is over 400 mg/dL (22.2 mmol/L) or below 40 mg/dL (2.2 mmol/L).
- **Manual:** A blood glucose value that was not entered as a calibration.

Exporting Data

From the Transmitter Log Report page, you can also export data in a .xlsx format.

1. Select Reports > Transmitter Log

2. Select “Export Data” at the top of each table to export the data as an .xlsx (Excel Spreadsheet)


Sensor Glucose					
Search					
Export Data					
Date	Time	Sensor Glucose	Unit	Sensor ID	
14-Feb-2025	06:55 PM	110	mg/dL	7679	
14-Feb-2025	06:50 PM	114	mg/dL	7679	
14-Feb-2025	06:45 PM	109	mg/dL	7679	
14-Feb-2025	06:40 PM	109	mg/dL	7679	

6. Settings

Glucose

Glucose settings allow you to customise the glucose target and alert levels in the reports. This setting impacts how glucose and other event data are organised and displayed in reports. When the patient wirelessly syncs their data via the mobile app, the glucose settings they have saved on the app will also be synced to their DMS account.


IMPORTANT: Glucose settings DO NOT get synced when data is uploaded from the smart transmitter via the USB Cable.



To change the settings in the Eversense DMS Pro view, follow the steps below:

1. On the top right of the Dashboard, click “Glucose Settings”
2. On the pop-up, click the mg/dL or mmol/L button below the Units Of Measure to select the desired unit of measure for displaying glucose data.
If you change the unit of measure, the glucose range and alert range values will change accordingly.
3. Click the +/- buttons to the right of the High Target and Low Target to display and adjust the high and low limits of your glucose Target Range
4. Repeat for the glucose Alert Range
5. When you are done making changes, click “Save” to save your changes

Note: These settings will be saved for the next time you view this patient’s reports, even after logging out of DMS Pro. Eversense DMS Pro settings do not affect the patients’ Eversense CGM System settings or the Eversense DMS display.




7. Help

About

This page shows the most up-to-date versions of the Eversense DMS Pro.

IMPORTANT: For questions about your Eversense DMS Pro account, contact your local distributor or visit global.eversenseddiabetes.com to locate your local distributor.

7



8. Troubleshooting

1. My previous filters are no longer applied upon going back to my Home page or Patients page.

- Filters are refreshed upon leaving the Home page or Patients tab. Please reapply filters to see the desired list of patients.

2. If the only clinic admin inadvertently removes their admin privileges please contact Customer Care.

3. CGM users in Germany may decline to create an Eversense account. These users cannot share data to Eversense DMS Pro.

Basic UDI/DI:

- Eversense DMS Pro: 081749102FG5800U5

Call 04-6175390 for Customer Support.

Distributed by:

DYN Diagnostics Ltd.

7 Ha'eshel St. P.O. Box 3063

Caesarea Industrial Park 3079504, Israel


Phone: 04-6175390

Email: patient.care@dyn.co.il



Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



 **Senseonics, Inc.**
20451 Seneca Meadows Parkway
Germantown, MD 20876-7005 USA
global.eversensedibetes.com

Patents:
www.senseonics.com/products/patents

